



# Baby Scale BBBYSCLXA Quick Start Guide

## Product Specifications

Work temperature: 10C ~ 40C  
Environmental Humidity: <80%  
Max. Capacity: 20kg  
Min. Capacity: 20g  
Graduation: 5g  
Power Supply: 4 x AA batteries  
Low Power: Lo  
Error Indicator: Err  
Measurement Units: kg or lb  
Memory Function: Optionally save and recall weight values  
Automatic Shut-off

## Keys and Functions

**ON/TARE:** Press to power the scale on and start weighing process. For tare operations, press the key a second time to clear the display and the “zero” annunciation will display in the bottom left corner of the display screen. If you press and hold this key for 3 seconds, the scale will turn off.

**MODE:** Weight unit toggle. Press this key to swap between KG and LB measurements in turn. The current weight unit will be shown on the display screen.

**MEM:** While scale is in a zero weight situation “0.00”, press this key to recall the last weight value saved in memory, and the corresponding indicator will illuminate. (If there is no previously saved data, the screen will show “Err”.) Press this key to save the current stable weight value while the current weight is displaying.

## Operations

Place the batteries in the battery compartment and press "ON/TARE" to start the scale. Wait until "0.00" is displayed.

Gently place the baby in the pan, and the weight will be displayed in the display screen after a few seconds. If the baby is left lying in the pan after its weight value has been shown, the weight will continue to display for 5 minutes before the unit turns itself off.

If the baby is removed from the pan after the weight has been displayed, the scale will return to "0.00" and will turn off automatically after 5 minutes.

To check the former saved value, press the "MEM" key while the scale is displaying "0.00". The scale will save the option value during each weight operation.

## Cautions

The scale should be placed on a steady, flat, hard surface to ensure the utmost safety for the baby.

You must replace the batteries immediately when the window displays "Lo".

Keep the scale clean to ensure normal working operations.

If the scale will not be used for an extended period of time, the batteries should be removed.

Please place the baby in the center of the pan to ensure the accuracy of the weight operation.

Do not place anything or anyone heavier than 20.5kg onto the scale.