

kogan

СОВЕТЫ ПО ЭКСПЛУАТАЦИИ

USER MANUAL



Thank you for buying the Kogan Deluxe Hand-held Massager, the ultimate aid in promoting healthy blood circulation and soothing aching muscles.



Product Specification

Description:	Deluxe Hand-held Massager
Model:	KAMSGHANDAA
Voltage:	220~240V
Frequency:	50~60Hz
Power:	25W
Safety Design:	<input type="checkbox"/>
Rated Working Time:	15 Minutes


Usage Positions



Name and usage of components



Split: Please press the split button on both sides, forward slightly harder then pull back.



Combination: Please insert the end lid into the tail of the massager and it can then be used in split mode. Pull the tail lid out and reassemble with the handle to be used in full massager mode.

Usage

- Please confirm the massager function switch is in the OFF position prior to inserting or removing the power plug from the wall socket.
- Put the function switch in (M) and adjust the Massage Speed Switch to suit your desired massage strength. Rotate clockwise and the massage strength will increase gradually. Rotate the Speed Switch counter-clockwise, and the massage strength will gradually decrease.
- Put the function switch into (M+L) state to get both massage and heating.
- It is not recommended that you massage the same region of the body for longer than 15 minutes at a time. Doing so can cause excess stimulation, which reduces the long term effect of the massage. We recommend that you massage for 10~15 minutes, then leave the region to rest before resuming further massage.
- When you have finished using the massager, turn the power OFF and then remove the plug from the power socket.

Caution

- Please ensure the machine is connected to the specific voltage listed in this manual.
- Should not be used on the torso of pregnant women.
- Do not massage over skin covering an infection or a bone fracture.
- This massager must not be used on a wet body.
- This massager must not be used by someone with a wet hand.
- Never keep this device in water or near the shower. Do not store the device in wet rooms such as the bathroom or laundry.
- Never use this device near combustible materials or heat sources.
- Never use this device during sleep.
- Never use this device on infants, toddlers or animals.
- Never store the device with the plug still attached to the power supply.
- Never use the massager for an extended period of time on the head, chest, knees or elbow joints as further pain will most likely result.
- Never allow children to use this massager unattended.
- Don't insert your hand or fingers into the gap between the massage head and the case, as this can result in pinching.
- Never carry the massager by the power cord.
- Never disassemble the device yourself, as this may cause a malfunction or electric shock. Any modifications performed by the user instantly voids the warranty and absolves Kogan of any legal liability for injuries or damage caused.
- Disconnect from the power supply when cleaning. Please use a neutral cleaner. Never use engine oil or organic solvents.