





Percussion Neck & Shoulder Massager KAMSGSHNECA

Product Layout

Massager



Control Panel



Operation

1. Connect the power plug to an electrical outlet and turn the power outlet on.
2. Place the massage unit onto the part of the body you want massaged, and make sure it is secured comfortably.
3. Press the POWER/MODE button to enter automatic massage mode.
4. Press the STRENGTH +/STRENGTH - buttons to adjust the intensity of the massage.
5. After 15 minutes the massager will stop functioning.
6. You can press the POWER/MODE button to cycle through the various massage modes.

Warnings

- It is recommended that you do not use the massager for longer than 15 minutes at a time.
- Stop use immediately if you feel any discomfort.
- Ensure the power plug has been disconnected from the power outlet when not in use.
- In the event of a power failure, disconnect the power cord to prevent sudden operation upon the power being restored.
- Do not let children play with the massager.

Attention

Please handle this product with care and inspect it regularly to ensure it is in good working order.

If the product, power supply cord or plug shows any signs of damage: stop use, unplug and contact Kogan.com support.