

ALL THE BENEFITS, NONE OF THE EFFORT

**HIGH-PERFORMANCE
SLIMMING UNIT**

**USER
MANUAL**

PG-2001G

- ▶ **Strong Vibration**
- ▶ **Oscillating & Screwing**
- ▶ **Detox & Massage**



Heat



Ventilating
System



Oscillating
Massager



5 Level

For use on a variety of areas



PG-2001G USER MANUAL

FEATURE

1. Oscillating for massage, two vibrate wheels designed to make slimming effective.
2. With two auto modes (automatically run a speed cycle from level 1 to level 5)
3. With a manual modes (speed and intensity available from level 1 to level 5 for selection)
4. With 5, 8, 10, 12 minutes auto off function.
5. Heating function

HOW TO OPERATE THE UNIT

1. Wrap the belt firmly to your body and adjust the size by the fastener.
2. Put the adaptor jack of the adaptor to the Power connector Jack of the belt.
3. Plug in the power adaptor to the wall socket. (Make sure that the input voltage of the wall socket) adaptor is suitable for the supplied adaptor.
4. To turn on the unit, press the on/off power button once, then the belt will be run the auto mode, or you can press the mode button to adjust to manual mode.
5. To choose auto manual, heat mode by pressing the function button with the green light is on.
6. To choose intensity by pressing "strong" or "weak" button is just valid under manual mode. Pressing button without any sounds represents the strongest mode.

7. To choose four kinds of auto timing: 5, 8, 10, 12 minutes with the red light is on. It auto stops. Press "on/off" to restart. Please do not use for more than 30 minutes one time.

PRECAUTION

1. Use whilst pregnant is not allowed.
2. Please stop using if abnormal situation happens.
3. When using for a long time, the product will stop waiting for the ventilating device making it cooler to a safe temperature.
4. Using the product, it's normal to feel itchy because of blood circulation speed up suddenly, tapping 1 to 3 minutes later, itchy feeling disappear.
5. When not in use, make sure to plug out the adaptor from the electricity socket.
6. Do not repair or maintain the product.
7. Do not use in hot and wet place.
8. Do not store the unit in a sunny place or a damp place for long time.
9. Do not put heavy articles on the unit.
10. For cleaning the belt, let the belt cool down and clean the dust with a dry cloth .

TECHNICAL SPECIFICATION

Model NO. PG-2001G

Name: High-performance Slimming Unit

Input voltage: DC 24V 1500mA

