

kügan

KATENELMASA

USER MANUAL



Features

The Nerve Simulation Massager offers the following features:

- ↓ Uses digital technology to ensure efficient massage.
- ↓ This unit is equipped with a large LCD that displays the intensity functions and massage time.
- ↓ Intensity buttons allow you to instantly change the magnitude of the output pulse.
- ↓ Numerous waveform options are available.

Main Unit



Accessories

The Nerve Simulation Massager comes with the following accessories:

- ↓ Electrode Pads
- ↓ 2-in-1 Electrode Wire

Specifications

Power	DC 4.5v, 3 AAA batteries
Frequency	1-330hz
Power Consumption	Less than 180mw
Battery Life	45 continuous hours under normal operating conditions
Device Size	153mm, 70mm, 18mm
Weight	90g

Instructions

Read the following instructions before using your Nerve Simulation Massager to ensure proper operation.

Getting Started

Install the Batteries

1. Remove the battery cover.
2. Insert three AAA batteries, paying attention to the polarity of the batteries.
3. Replace the cover and close firmly.

Prepare the Electrodes

1. Attach the adhesive pad to the electrode cord.
2. Remove the protective film from the adhesive pads shortly prior to using the device.

Place the Electrodes

1. Clean the portion of the body that will be massaged with a damp towel to remove oil, dirt, sweat, and make up.
2. Gently stick the electrodes directly onto the skin that will be massaged. (**Note:** Massage will not be effective without two electrodes placed on the body; however, it is not necessary to place the two electrodes symmetrically.)

Using the Nerve Simulation Massager

Complete the following steps to begin using your Nerve Simulation Massager.

1. Press the **On/Off** button to turn on the unit. The LCD screen will display the power status, as well as the current mode, intensity, and duration. The default mode is acupuncture, with an intensity of 0 and a duration of 15 minutes.
2. Press the **En/Ch** button to select the language (English or Chinese) for the display.
3. Press the **Mode** button to select from the eight available modes (acupuncture, stroke, massage, cupping, manipulation, scraping, weight reducing, immunotherapy).
4. Press the **Strong** button to increase the intensity level of the massage by one level or press the **Weak** button to decrease the intensity level of the massage by one level. Fifteen intensity levels are available.
5. Press the **Time** button to select the length of the massage (in increments of five minutes). Fifteen minutes is the maximum duration for a single massage.

Safety Precautions

Equipment damage or danger to the user can occur if you use the Nerve Simulation Massager in conjunction with electronic medical equipment. Do **NOT** use the Nerve Simulation Massager if you use any of the following:

- ⌋ A heart rate regulator or any other internal electronic medical device.
- ⌋ An artificial cardiopneumatic pacemaker, electronic life support, or other internal medical device.
- ⌋ A cardiometer or other wearable electronic medical device.

Speak to your doctor before using the Nerve Simulation Massager if you are under medical treatment or have any of the following conditions:

- ⌋ Malignant tumour
- ⌋ VD infection
- ⌋ Menstruating or pregnant
- ⌋ Heart disease
- ⌋ Fever above 39°
- ⌋ Abnormal cutaneous anaphylaxis
- ⌋ Area of disease near your eyes, inside your mouth, on the genitalia or near points of external skin trauma.

Do **NOT** use the Nerve Simulation Massager simultaneously with other instruments or in conjunction with lubricants, lotions or oils.

For safety purposes, if you change the mode during a massage session, the intensity will automatically decrease to 0. You can readjust the intensity to the appropriate level after the new mode is set.

If you want to change the position of an electrode during a massage session, turn off the device first. You can restart after you have repositioned the electrode.

When removing the electrodes from the device, always pull from the body of the device, not the cord.